

FIG. 1




FIG. 2




FIG. 3

41

SCOPE TYPE NAME	SCOPE ID	USE START DATE	NUMBER OF USE TIMES	USE TIME

FIG. 4A

42

USE DATE	USE START TIME	USER	PATIENT ID	PATIENT NAME	EXAMINATION ITEM	CLEANED/NOT CLEANED

FIG. 4B




FIG. 5





FIG. 6



**FIG. 7**

STATISTICS/HISTORY

DAILY REPORT

MONTHLY REPORT

PERFORMANCE LIST

OUTPUT

TABLE FILE

OUTPUT

CLEANING HISTORY

SCOPE HISTORY

FOLLOW-UP SEARCH

SCOPE LIST

UPDATE

CLASSIFICATION: ALL ARE DISPLAYED

PRINT LIST

SCOPE NO	CLASSIFICATION	SCOPE NAME	SCOPE ID	USE START DATE	NUMBER OF USE TIMES	USE TIME (MINUTES)	CURRENT STATE	LIFE
01	GIF	GIF-H260	100000001	2001/08/20	10	125	CLEANED	<div></div>
02	GIF	GIF-GX240	100000002	2001/08/20	30	300	CLEANED	<div></div>
03	GIF	GIF-KQ240	100000003	2001/08/20	15	154	CLEANED	<div></div>
04	GIF	GIF-XK240	100000004	2001/08/20	3	37	CLEANED	<div></div>
05	CF	CF-Q240A1	200000001	2001/08/20	5	100	CLEANED	<div></div>
06	CF	CF-24001	200000002	2001/08/20	10	120	NOT CLEANED	<div></div>
07	EUS	CF-UNQ240	300000001	2001/08/20	1	15	NOT CLEANED	<div></div>
08	ERCP	TJF-240	400000001	2001/08/20	3	45	NOT CLEANED	<div></div>
09	BF	BF-240	500000001	2001/08/20	2	53	NOT CLEANED	<div></div>

LOGOUT

LOGIN USER:

FIG. 8



# SCOPE USE HISTORY

■ SCOPE  
INFORMATION

SCOPE TYPE  
NAME :  
GIF-H260

SCOPE ID :  
10000001

USE START DATE :  
2002/08/28

NUMBER OF USE  
TIMES :  
10 TIMES

USE TIME  
(MINUTES) :  
125 MINUTES

SCOPE LIFE



■ USE HISTORY

DISPLAY DURATION : 2002/08/29 ~ 2002/08/30 USER : ALL

USE DATE	USE START TIME	USER	PATIENT ID	PATIENT NAME	EXAMINATION ITEM	CLEANED/ NOT CLEANED
2002.08.29	09:00:00	OLYMPUS TARO	00140481	OLYMPUS ICIRO	COLONOSCOPY	CLEANED
2002.08.29	09:30:00	OLYMPUS TARO	00000001	OLYMPUS JIRO	COLONOSCOPY	CLEANED
2002.08.30	11:35:30	OLYMPUS TARO	11112222	OLYMPUS SABURO	ERCP	CLEANED
2002.08.30	14:10:00	OLYMPUS TARO	99999999	OLYMPUS SHIRO	UPPER PORTION EMR	CLEANED

BACK TO LIST

LOGIN USER

FIG. 9




FIG. 10